



Wal-Mart Inspection Form

Today's Date: _____

Company Name: Wal-Mart	Store #	Address/ City/ State:	
Inspectors Name:		AED Cabinet Physical Location (Be Specific):	
Manager's Name:	AED Model: CR2	Serial Number:	

Check the following	PASS	FAIL	Comments
Is the Readiness Indicator Green or ready to use?			
Is unit reporting connected to Wi-Fi during two-button test?			
Are electrodes (Pads) connected? Verify expiration date to the right			Exp. Date:
Are Batteries connected properly? Verify expiration date to the right			Exp. Date:
Is the unit clean, undamaged, and free of excessive wear / cracks?			
When lid is closed, the unit should NOT beep, chirp, or continue talking.			
Is response kit attached or stored next to unit inside cabinet?			
Cabinet and wall sign were installed according to directions?			
Cabinet alarm operates correctly?			

Comments:

UPLOAD FORM VIA THE CARDIO PARTNERS PORTAL
 ANY WIFI ISSUES, QUESTIONS, OR CONCERNS PLEASE CALL
 866-349-4363 ext. 6
 or
 Christian Hernandez @ 214-501-5104